

Psychopy: a Python library for psychological experiments.

Session 2 – PsychoPy Introduction

Jaime A. Riascos Salas

**Computational Neurosciences and Cybernetics
Systems (CONCYS)**

Institute of Informatics

Federal University Rio Grande do Sul (UFRGS)

Porto Alegre – Brazil.

- 1. Introduction PsychoPy**
- 2. Installation**
- 3. Architecture – Builder and Coder**
- 4. Understanding PsychoPy**
- 5. Monitor object**
- 6. Stimuli objects**
- 7. Hands-on!!!**



PsychoPy

What is it ?

It is an open-source package for creating and running psychological experiments. PsychoPy uses both OpenGL (most widely adopted 2D and 3D graphics library) and Python to give scientists and students a free and simple stimulus presentation and control package. It is widely used by many labs worldwide for psychophysics, cognitive neuroscience and experimental psychology.

www.psychopy.org



PsychoPy

Advantages?

Instead the paid software (Matlab, up 600 \$USD; e-Prime, up 300 \$USD), PsychoPy is totally free and open-source.

As any open source software, you are able to download and modify the package as much as you want. PsychoPy is developed in a single piece of software enabling to be precise enough for psychophysics, intuitive enough for undergraduate psychology and flexible enough for anyone who is intending to develop own experiments.



PsychoPy

Features?

- Simple install process;
- Precise timing;
- Huge variety of stimuli (see screenshots) generated in real-time:
 - Images
 - Random dots
 - Text (unicode in any truetype font)
 - Shapes and movies (DivX, mov, mpg...)
 - sounds (tones, numpy arrays, wav, ogg...)



PsychoPy

Features?

- Platform independent - run the same script on Win, macOS or Linux;
- Flexible stimulus units (degrees, cm, or pixels);
- Input from keyboard, mouse, microphone or button boxes;
- Multi-monitor support;
- Automated monitor calibration (for supported photometers);
- **Coder interface**, for those that like to program
- **Builder interface**, for those that don't



PsychoPy - Installation

Windows/Mac users

Linux users

Install the Standalone package:

Code Issues 166 Pull requests 4 Projects 0 Wiki Insights

Releases Tags

Pre-release

3.0.0b3
1af81bf

peircej released this 4 hours ago · 0 commits to master since this release

Release 3.0.0b3

Assets

PsychoPy-3.0.0b3.zip	11 MB
StandalonePsychoPy3-3.0.0b3-MacOS.dmg	300 MB
StandalonePsychoPy3-3.0.0b3-win32.exe	204 MB
StandalonePsychoPy3_PY3-3.0.0b3-MacOS.dmg	328 MB
StandalonePsychoPy3_PY3-3.0.0b3-win32.exe	244 MB
Source code (zip)	
Source code (tar.gz)	

- `sudo apt-get install psychopy`
(maybe you need to install some dependencies)

- Neurodebian
(<http://neuro.debian.net/>)



PsychoPy - Installation

Manual Installation (win, linux, mac)

Now that most python libraries can be install using `pip` it's relatively easy to manually install PsychoPy and all it's dependencies to your own installation of Python. That isn't the officially-supported method (because we can't track which versions of packages you have) but for many people it's the preferred option if they use Python for other things as well.

Dependencies

You need a copy of Python 2.7.x from here, wxPython and probably pyo (or use an alternative audio library listed below). None of these support `pip install` yet so you need to download them:

- Python itself: <http://www.python.org/download/> (**version 3.x is not supported yet**)
- wxPython: <https://wxpython.org/download.php>
- pyo audio: <http://ajaxsoundstudio.com/software/pyo/>
- PyQt4 or PyQt5 are handy but not required and need manual installation

Then, if you want **everything** available you could paste this in to your terminal/commandline and go and get a coffee (will take maybe 20mins to download and install everything?):

```
pip install numpy scipy matplotlib pandas pyopengl pyglet pillow moviepy lxml openpyxl xlrd configobj  
pip install pyserial pyparallel egi iolabs  
pip install pytest coverage sphinx
```

Needed on Windows:

```
pip install pypiwin32
```

Needed on macOS:

```
pip install pyobjc-core pyobjc-framework-Quartz
```

<http://psychopy.org/installation.html>



PsychoPy Builder

For non-programmer users

PsychoPy Coder

For programmer users

PsychoPy Library

Stimuli

Timing

Trial

**Packing like a
Matlab Toolbox**

Python

**Programming
interface**

OpenGL


**Maths
(Numpy)**

**User
Interface
(wxPython)**

Low-level libraries



PsychoPy Builder





PsychoPy Builder

Components


A right-placed panel where you can find all objects for creating stimulus and adding an input/output device. A default list of components content:

- Aperture Component
- Cedrus Button Box Component
- Code Component
- Dots (RDK) Component
- Grating Component
- Image Component
- ioLab Systems buttonbox Component
- Keyboard Component
- Text Component
- Variable Component
- Microphone Component
- Mouse Component
- Movie Component
- Parallel Port Out Component
- Patch (image) Component
- Polygon (shape) Component
- RatingScale Component
- Sound Component
- Static Component





PsychoPy Builder




A middle-placed panel where you can find any number of Routines used in the experiment, describing the timing of stimuli, instructions and responses. These are portrayed in a simple track-based view, similar to that of video-editing software, which allows stimuli to come on go off repeatedly and to overlap with each other.




PsychoPy Builder

Flow

A bottom-placed panel where you can control all routines and loops used to form an experiment. All experiments have exactly one Flow. This takes the form of a standard flowchart allowing a sequence of routines to occur one after another, and for loops to be inserted around one or more of the Routines. The loop also controls variables that change between repetitions, such as stimulus attributes.




PsychoPy Builder





What are the main steps to create any experiment using PsychoPy?

PsychoPy uses a monitor for presenting stimulus (text, images, sounds, videos); thus, how these objects are showed on the monitor (duration and way) is known as routines, therefore, all of the created routines represent the whole experiment itself. Within the routines, how the data (user answers) is collected should be set-up in order to save this information for a posterior analysis. 1234





Monitor Center

PsychoPy provides a simple and intuitive way for you to calibrate your monitor and provide other information about it and then import that information into your experiment.

Information is inserted in the Monitor Center (Tools menu), which allows you to store information about multiple monitors and keep track of multiple calibrations for the same monitor.

Also some attributes about the monitor can be set on Experiment Settings (next slides)





Monitor Center



PsychoPy Coder

PsychoPy Builder

```
from psychopy import visual
win = visual.Window([1024,768], mon='SonyG500')
```

Class **psychopy.visual.Window**(size=(800, 600), pos=None, color=(0, 0, 0), colorSpace='rgb', rgb=None, dkl=None, lms=None, fullscr=None, allowGUI=None, monitor=None, bitsMode=None, winType=None, units=None, gamma=None, blendMode='avg', screen=0, viewScale=None, viewPos=None, viewOri=0.0, waitBlanking=True, allowStencil=False, multiSample=False, numSamples=2, stereo=False, name='window1', checkTiming=True, useFBO=False, useRetina=True, autoLog=True, *args, **kwargs)

PsychoPy2 Monitor Center

File Edit

Choose Monitor

testMonitor

New...

Save

Delete

Copy...

Delete

2018_07_23 22:08

Monitor Info

☐ Use Bits++

Screen Distance (cm): 57

Size (pixels; Horiz,Vert): 1024 768

Screen Width (cm): 30

Calibration Date: 2018 07 23 22:08

Notes:

default (not very useful) monitor

Calibration

Get Photometer

Scan all ports

Gamma Calibration...

Gamma Test...

Plot gamma

Chromatic Calibration...

Plot spectra

Linearization

easy: $a+kx^g$

	Min	Max	Gamma	a	b	k
lum	0.0000	1.0000	1.0000	nan	nan	nan
R	0.0000	1.0000	1.0000	nan	nan	nan
G	0.0000	1.0000	1.0000	nan	nan	nan
B	0.0000	1.0000	1.0000	nan	nan	nan

LMS->RGB

	L	M	S
R	0.0000	0.0000	0.0000
G	0.0000	0.0000	0.0000
B	0.0000	0.0000	0.0000

DKL->RGB

	Lum	L-M	L+M-S
R	0.0000	0.0000	0.0000
G	0.0000	0.0000	0.0000
B	0.0000	0.0000	0.0000



Experiment settings



PsychoPy Coder

PsychoPy Builder

```
from psychopy import visual
win = visual.Window([1024,768], mon='SonyG500')
```

Class **psychopy.visual.Window**(size=(800, 600), pos=None, color=(0, 0, 0), colorSpace='rgb', rgb=None, dkl=None, lms=None, fullscr=None, allowGUI=None, monitor=None, bitsMode=None, winType=None, units=None, gamma=None, blendMode='avg', screen=0, viewScale=None, viewPos=None, viewOri=0.0, waitBlanking=True, allowStencil=False, multiSample=False, numSamples=2, stereo=False, name='window1', checkTiming=True, useFBO=False, useRetina=True, autoLog=True, *args, **kwargs)

The screenshot shows the 'Experiment Settings' dialog box with the 'Screen' tab selected. The settings are as follows:


- Monitor: testMonitor
- Screen: 1
- Full-screen window: ☐
- Window size (pixels): [400,400]
- Color: \$[-0.004,-0.004,-0.004]
- Color space: rgb
- Units: use preferences
- Show mouse: ☐
- Blend mode: average

Buttons at the bottom: Help, OK, Cancel.



Stimulus


- Aperture Component
- Cedrus Button Box Component
- Code Component
- Dots (RDK) Component
- Grating Component
- Image Component
- ioLab Systems buttonbox Component
- Keyboard Component
- Microphone Component
- Mouse Component
- Movie Component
- Parallel Port Out Component
- Patch (image) Component
- Polygon (shape) Component
- RatingScale Component
- Sound Component
- Static Component
- Text Component
- Variable Component





Stimulus

- Aperture Component
- Cedrus Button Box Component
- Code Component
- Dots (RDK) Component
- Grating Component
- **Image Component**
- ioLab Systems buttonbox Component
- Keyboard Component
- Microphone Component
- Mouse Component
- Movie Component
- Parallel Port Out Component
- Patch (image) Component
- Polygon (shape) Component
- RatingScale Component
- **Sound Component**
- Static Component
- **Text Component**
- Variable Component





Stimulus

- Aperture Component
- Cedrus Button Box Component
- Code Component
- Dots (RDK) Component
- Grating Component
- **Image Component**
- ioLab Systems buttonbox Component
- Keyboard Component
- Microphone Component
- Mouse Component
- Movie Component
- Parallel Port Out Component
- Patch (image) Component
- Polygon (shape) Component
- RatingScale Component
- **Sound Component**
- Static Component
- **Text Component**
- Variable Component





Setting stimulus attributes

Stimulus attributes are features that each object has, for example, in a text, we have color, position, size... (Later)

Typically, they are set using either:

- A string, which is just some characters (a `message.setText('world')`)
- A scalar (a number)
- An x,y-pair (two numbers)



Setting stimulus attributes

Real world units

One of the particular features of PsychoPy is that you can specify the size and location of stimuli in units that are independent of your particular setup, such as degrees of visual. In order for this to be possible you need to inform PsychoPy of some characteristics of your monitor. Your choice of units determines the information you need to provide:

Units	Requires
'norm' (normalised to width/height)	n/a
'pix' (pixels)	Screen width in pixels
'cm' (centimeters on the screen)	Screen width in pixels and screen width in cm
'deg' (degrees of visual angle)	Screen width (pixels), screen width (cm) and distance (cm)



Text stimuli

PsychoPy Coder

PsychoPy Builder



```
from psychopy import visual

message = visual.TextStim(win, text='hello')
```

```
class psychopy.visual.TextStim(win, text='Hello
World', font="", pos=(0.0, 0.0), depth=0, rgb=None,
color=(1.0, 1.0, 1.0), colorSpace='rgb', opacity=1.0,
contrast=1.0, units="", ori=0.0, height=None,
antialias=True, bold=False, italic=False,
alignHoriz='center', alignVert='center', fontFiles=(),
wrapWidth=None, flipHoriz=False, flipVert=False,
languageStyle='LTR', name=None, autoLog=None)
```

text Properties

Basic Advanced

Name

Start Expected start (s)

Stop Expected duration (s)

Color

Font

Letter height \$

Position [x,y] \$

Text

Help OK Cancel



Image stimuli

PsychoPy Coder

PsychoPy Builder

```
from psychopy import visual

img1 = visual.ImageStim(win=win,
                        image=filename,
                        units='pix',
                        size=(200,200))
```

class **psychopy.visual.ImageStim**(win, image=None, mask=None, units="", pos=(0.0, 0.0), size=None, ori=0.0, color=(1.0, 1.0, 1.0), colorSpace='rgb', contrast=1.0, opacity=1.0, depth=0, interpolate=False, flipHoriz=False, flipVert=False, texRes=128, name=None, autoLog=None, maskParams=None)



<http://www.psychopy.org/api/visual/imagestim.html#psychopy.visual.ImageStim>



Audio stimuli

PsychoPy Coder

PsychoPy Builder



```
from psychopy import sound

aud1 = sound.Sound(value=soundfile,
                    volume=1)
```

```
class
psychopy.sound.backend_sounddevice.S
oundDeviceSound(value='C', secs=0.5,
octave=4, stereo=-1, volume=1.0, loops=0,
sampleRate=44100, blockSize=128,
preBuffer=-1, hamming=True, startTime=0,
stopTime=-1, name="", autoLog=True)
```



sound_1 Properties

Basic

Name


Start 0.0

Stop 1.0

Sound

Volume

Our first PsychoPy code!!!



**PsychoPy
Builder**

**Let's do it our first window with a stimuli
(text) using the builder.**

- Please open PsychoPy (using the Virtual Machine);
- Set up the window to be used with a size of 400x400 and black as background color;
- Create a stimuli text with the word “Hello”, color white and duration of two seconds;
- Create another stimuli text with the word “World”, the same color and duration of two seconds. This must start after the previous one.

Lets go into hacking!!!



PsychoPy
Coder

Let's do the previous example but this time using Python coding (Coder)

- Please open Spyder or any editor text (using the Virtual Machine);
- Call the necessary libraries (PsychoPy)
- Set up the window to be used with a size of 400x400 and black as background color;
- Create a stimuli text object with the word "Hello", color white;
- Set the draw operation;
- Call flip function (This really do the changes);
- Set the duration;
- Set the new value of the text object ("World");
- Repeat the same instructions after draw operation;

Ok, Cool, but only text is pretty bored!

Use some images!!!



PsychoPy
Builder

Let's do the previous example but this time using two images (Coder)

- Please open PsychoPy (using the Virtual Machine);
- Set up the window to be used with a size of 400x400 and white as background color;
- Create a stimuli image which will charge the file “hello.png”, with a duration of two seconds;
- Create another stimuli image which will charge the file “world.png” and duration of two seconds. This must start after the previous one.

The files can be found in the /media folder.

Ok Drake! Let's coding this example



PsychoPy
Builder



PsychoPy
Coder

- Please open Spyder or any editor text (using the Virtual Machine);
- Call the necessary libraries (PsychoPy)
- Set up the window to be used with a size of 400x400 and white as background color;
- Create an image object which has the file “hello.png”;
- Create another image object which has the file “world.png”;
- Repeat the same procedure as in the text stimuli.

The files can be found in the /media folder.

But... so for what do we show an image if actually, we can use an audio?

PsychoPy
Builder



NOT BAD

- Please open Spyder or any editor text (using the Virtual Machine);
- Call the necessary libraries (PsychoPy)
- Set up the window to be used with a size of 400x400 and gray as background color;
- Create an audio object which has the file “hello.wav” with a duration of one second;
- Create another audio object which has the file “world.wav” with a duration of one;

The files can be found in the /media folder.



Coding again!!

- Please open Spyder or any editor text (using the Virtual Machine);
- Call the necessary libraries (PsychoPy)
- Set up the window to be used with a size of 400x400 and gray as background color;
- Create an audio object which has the file “hello.wav”;
- Create another audio object which has the file “world.wav”;
- Repeat the same procedure as in the image stimuli.

The files can be found in the /media folder.

PsychoPy
Coder

PsychoPy
Builder

Let's do the previous examples all together!!!



- Please open both Builder and Spyder or any editor text (using the Virtual Machine);
- Set up the window to be used with a size of 800x800 and white as background color;
- Create an audio, image and text object with the “Hello” instance.
- The image must be placed at the centre and the text upper.
- All objects should appear together.
- The “world” objects (image, text and audio) should appear after.

[1] <http://psychopy.org/documentation.html>

**Any
questions... ?**

